





HERBERT
A new Sense of Confidence

PECTUS BAR SET



GW2ZXBT225M

Herbert (Thailand) Co., Ltd.

47 Moo 1, Samran Sub-district, Mueang District, Khon Kaen 40000, Thailand

Website : www.herbert.co.th Email : Info@herbert.co.th



Pectus Bar Set



Pectus Tunneling - Short
GW2ZXBT225D



Pectus Tunneling - long
GW2ZXBT225E



Pectus Handle Bender
GW2ZXBT225F



Bending Cylinder 9.7 mm
GW2ZXBT225G



Bending Cylinder 11.7 mm
GW2ZXBT225H



Bending Cylinder 13.7 mm
GW2ZXBT225J



Pectus Lever Bender
GW2ZXBT225K



Pectus Table Top Bender
GW2ZXBT225L




Pectus Box
GW2ZXBT2221

Specifications

	Code	Dimension	Material
Pectus Tunneling - Short	GW2ZXBT225D	20 x 451.3 x 81.6 mm	-
Pectus Tunneling - Long	GW2ZXBT225E	20 x 484.3 x 88.6 mm	-
Pectus Handle Bender	GW2ZXBT225F	118 x 205 x 29 mm	Stainless Steel SUS420
Bending Cylinder 9.7 mm	GW2ZXBT225G	Ø9.7 x 19.5 mm	Stainless Steel SUS420
Bending Cylinder 11.7 mm	GW2ZXBT225H	Ø11.7 x 19.5 mm	Stainless Steel SUS420
Bending Cylinder 13.8 mm	GW2ZXBT225J	Ø13.8 x 19.5 mm	Stainless Steel SUS420
Pectus Lever Bender	GW2ZXBT225K	36x267.3 x 56.5 mm	-
Pectus Table Top Bender	GW2ZXBT225L	80 x 410 x 340 mm	-
Pectus Box	GW2ZXBT2221	337 x 567 x 101 mm	Aluminium 5083


PECTUS BAR STABILIZER

Sliding Pectus Bar Stabilizers onto the Pectus Bar prevents the implanted Pectus Bar from tilting.



Pectus bar stabilizer, Low profile

GD253DMZ02



Specifications

	Code	Dimension	Material
Pectus bar stabilizer, Low profile	GD253DMZ02	15 x 60 x 5.3 mm	Stainless Steel 316 LVM

PECTUS BAR TEMPLATE

- Surgical preparation outside the body for determining a suitable implant
- Pectus Bar Templates may not be implanted




Specifications

	Code	Dimension	Material
Pectus Bar Template 178 mm	GW2ZXBT2252	12.7 x 178 x 1 mm	Aluminium 1100
Pectus Bar Template 203 mm	GW2ZXBT2253	12.7 x 203 x 1 mm	Aluminium 1100
Pectus Bar Template 229 mm	GW2ZXBT2254	12.7 x 229 x 1 mm	Aluminium 1100
Pectus Bar Template 254 mm	GW2ZXBT2255	12.7 x 254 x 1 mm	Aluminium 1100
Pectus Bar Template 279 mm	GW2ZXBT2256	12.7 x 279 x 1 mm	Aluminium 1100
Pectus Bar Template 305 mm	GW2ZXBT2257	12.7 x 305 x 1 mm	Aluminium 1100
Pectus Bar Template 330 mm	GW2ZXBT2258	12.7 x 330 x 1 mm	Aluminium 1100
Pectus Bar Template 356 mm	GW2ZXBT2259	12.7 x 356 x 1 mm	Aluminium 1100
Pectus Bar Template 381 mm	GW2ZXBT225A	12.7 x 381 x 1 mm	Aluminium 1100
Pectus Bar Template 406 mm	GW2ZXBT225B	12.7 x 406 x 1 mm	Aluminium 1100
Pectus Bar Template 432 mm	GW2ZXBT225C	12.7 x 432 x 1 mm	Aluminium 1100

PECTUS BAR

Serrated ends for secure attachment of a Pectus Bar Stabilizers
 Ø 5 mm drill holes on both sides for insertion of a thread-through band.



Specifications

	Code	Dimension	Material
Pectus Bar 203 mm	GD253DMZ031	12.7 x 203 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 216 mm	GD253DMZ032	12.7 x 216 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 229 mm	GD253DMZ033	12.7 x 229 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 241 mm	GD253DMZ028	12.7 x 241 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 254 mm	GD253DMZ029	12.7 x 254 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 267 mm	GD253DMZ024	12.7 x 267 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 279 mm	GD253DMZ018	12.7 x 279 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 293 mm	GD253DMZ019	12.7 x 293 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 305 mm	GD253DMZ020	12.7 x 305 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 318 mm	GD253DMZ021	12.7 x 318 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 330 mm	GD253DMZ022	12.7 x 330 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 343 mm	GD253DMZ025	12.7 x 343 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 356 mm	GD253DMZ026	12.7 x 356 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 368 mm	GD253DMZ030	12.7 x 368 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 381 mm	GD253DMZ027	12.7 x 381 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 394 mm	GD253DMZ034	12.7 x 394 x 3 mm	Stainless Steel 316 LVM
Pectus Bar 406 mm	GD253DMZ035	12.7 x 406 x 3 mm	Stainless Steel 316 LVM